

Degradome protocol

By ZZ & Bo, modified by Xuan 02-23-16, modified by Ildar 2021, modified by Ayca 2024

5' adapter (equimolar mix of BA5-UMI-a and BA5-UMI-b, 41-nt custom RNA oligos, PAGE purified):

5'-GUUCAGAGUUCUACAGUCCGACGAUCNNNCGANNNUACNNN-3'

5'-GUUCAGAGUUCUACAGUCCGACGAUCNNNAUCNNNAGUNNN-3'

DEG-RT (same as CAGE RT primer)

5'-GCACCCGAGAATTCCANNNNNNNN-3'

DEG-PCR-1L (same as CAGE PCR primer)

5'-CTACACGTTTCAGAGTTCTACAGTCCGA-3'

DEG-PCR-1R

5'-GCCTTGGCACCCGAGAATTCCA-3'

PCR primer1/BPCR1 (same as small RNA-Seq library PCR primer)

5'-AATGATACGGCGACACCGAGATCTACACGTTTCAGAGTTCTACAGTCCGA-3'

PCR primer2/BPCRIdX:

PCRId1: 5'-CAAGCAGAAGACGGCATACGAGATCGTGATGTGACTGGAGTTCCTTGGCACCCGAGAATTCCA-3'

PCRId2: 5'-CAAGCAGAAGACGGCATACGAGATACATCGGTGACTGGAGTTCCTTGGCACCCGAGAATTCCA-3'

PCRId3: 5'-CAAGCAGAAGACGGCATACGAGATGCCTAAGTGACTGGAGTTCCTTGGCACCCGAGAATTCCA-3'

PCRId4: 5'-CAAGCAGAAGACGGCATACGAGATTGGTCAAGTGACTGGAGTTCCTTGGCACCCGAGAATTCCA-3'

PCRId5: 5'-CAAGCAGAAGACGGCATACGAGATCACTGTGTGACTGGAGTTCCTTGGCACCCGAGAATTCCA-3'

PCRId6: 5'-CAAGCAGAAGACGGCATACGAGATATTGGCGTGACTGGAGTTCCTTGGCACCCGAGAATTCCA-3'

PCRId7: 5'-CAAGCAGAAGACGGCATACGAGATGATCTGGTGACTGGAGTTCCTTGGCACCCGAGAATTCCA-3'

PCRId8: 5'-CAAGCAGAAGACGGCATACGAGATTCAAGTGTGACTGGAGTTCCTTGGCACCCGAGAATTCCA-3'

PCRId9: 5'-CAAGCAGAAGACGGCATACGAGATCTGATCGTGACTGGAGTTCCTTGGCACCCGAGAATTCCA-3'

PCRId10: 5'-CAAGCAGAAGACGGCATACGAGATAAGCTAGTGACTGGAGTTCCTTGGCACCCGAGAATTCCA-3'

PCRId11: 5'-CAAGCAGAAGACGGCATACGAGATGTAGCCGTGACTGGAGTTCCTTGGCACCCGAGAATTCCA-3'

PCRId12: 5'-CAAGCAGAAGACGGCATACGAGATTACAAGGTGACTGGAGTTCCTTGGCACCCGAGAATTCCA-3'

PCRId13: 5'-CAAGCAGAAGACGGCATACGAGATTTGACTGTGACTGGAGTTCCTTGGCACCCGAGAATTCCA-3'

PCRId14: 5'-CAAGCAGAAGACGGCATACGAGATTTGACTGTGACTGGAGTTCCTTGGCACCCGAGAATTCCA-3'

PCRId15: 5'-CAAGCAGAAGACGGCATACGAGATTGACATGTGACTGGAGTTCCTTGGCACCCGAGAATTCCA-3'

PCRId16: 5'-CAAGCAGAAGACGGCATACGAGATGGACGGGTGACTGGAGTTCCTTGGCACCCGAGAATTCCA-3'

PCRId17: 5'-CAAGCAGAAGACGGCATACGAGATCTCTACGTGACTGGAGTTCCTTGGCACCCGAGAATTCCA-3'

PCRId18: 5'-CAAGCAGAAGACGGCATACGAGATGCGGACGTGACTGGAGTTCCTTGGCACCCGAGAATTCCA-3'

PROCEDURE

I. Ribosomal RNA depletion (also see rRNA depletion folder in Common Folder) & DNase treatment

Mouse: Pooled 186 rRNA antisense oligos (100 μ M/each) that tile 18S, 28S, 5S, 5.8S, mt 12S, and mt 16S rRNAs. Final working concentration: 0.5 μ M for each oligo.

Fly: Pooled 160 rRNA antisense oligos (100 μ M/each) that tile 18S, 28S, 5.8S, mt L, and mt S rRNAs. Final working concentration: 0.5 μ M for each oligo.

- Use 500 ng – 2 μ g high quality total RNA (strongly recommended) as starting material. Suggested kit for RNA extraction: mirVana by Ambion.
- Add 2 μ L of pooled rRNA AS oligos for every 1 μ g of total RNA (ratio: 1 μ L oligos for 1 μ g RNA).
- Make final volume 14 μ L with 1X rRNA oligo hybridization buffer (aliquots in the -20°C radioactivity common freezer).
- Heat the mixture to 95°C , then slowly cool it down ($-0.1^{\circ}\text{C}/\text{sec}$) to 22°C in a PCR machine. Incubate at 22°C for additional 5 minutes. Place on ice.
- Add the following to RNA (final volume: 20 μ L):
 - 2 μ L 10X RNase H digestion buffer
 - 2 μ L water
 - 2 μ L Thermostable RNase H (10 units; Hybridase from Epicentre)
- Incubate at 45°C for 30 minutes. Place on ice.
- *DNase treatment:* Add the following to each tube (final volume 50 μ L):
 - 23 μ L water
 - 5 μ L Turbo DNase 10X buffer
 - 2 μ L Turbo DNase
- Then incubate at 37°C for 20 minutes.
- Purify RNA using RNA Clean & Concentrator-5 (Zymo Research: Cat R1015), enriching RNA > 200nt and removing 5S rRNA and tRNAs. Note that table top centrifuge can be used for steps 3,4,5 and 7.
 1. Adjust *RNA Binding Buffer*: combine 50 μ L RNA binding buffer and 50 μ L 100% ethanol
 2. Add 100 μ L of adjusted *RNA Binding Buffer* to RNA samples. Mix well.
 3. Transfer the mixture to the Zymo-Spin column. Spin @ 17,000Xg for 30 seconds (RNA > 200nt retained in the column). Discard flow-through.
 4. Add 400 μ L *RNA Prep Buffer*. Spin @ 17,000Xg for 30 seconds. Discard flow-through.
 5. Add 700 μ L *RNA Wash Buffer*. Spin @ 17,000Xg for 30 seconds. Discard flow-through.
 6. Add 400 μ L *RNA Wash Buffer*. Spin @ 17,000Xg for 2 minutes. Discard flow-through
 7. Spin @ 17,000Xg for 30 seconds to remove any ethanol residue
 8. Move the column to a clean 1.7mL tube.
 9. Add 7 μ L water to the column. Incubate @ room temp for 1 minute. Spin @ 17,000Xg for 30 seconds to elute RNA. Transfer 6.5 μ L to a new tube.

II. Library preparation

Note: This protocol is written for using 2 μ g mouse testis/Hi5 total RNA as a starting point. If you work with a different organism, tissue, and/or amount of starting total RNA, it's important to optimize the PCR cycle numbers for both PCR steps.

1. 5' adapter ligation

- Add to the 6.5 μ L RNA from the previous step:

10x 5' RNA ligase buffer	2ul
50% PEG (20% final)	8ul
25uM BA5-UMI adapter	1ul
RNasein(40U/ul) (N2615)	0.5ul
T4 RNA Ligase High Conc(M0437M)	2ul

- Incubate at 16°C for 16 hours.
- *Optional:* Add 30 μ L water to the ligation reaction for a larger volume that's easier to handle.
- Purify RNA with RNA Clean & Concentrator-5 (Zymo Research: Cat R1015), which enriches for $\geq 200\text{nt}$ RNAs. Elute with 11 μ L water. Transfer 10 μ L RNA to a 200 μ L PCR tube.

2. First strand cDNA synthesis

- Add to the RNA from the previous step:

5X FS buffer	5 μ L
Deg-RT primer (50 μ M)	1 μ L

- 65°C for 3 min for 4°C forever. Place on ice.

- Add the following 9 μ L mixture:

Water	5.25 μ L
dNTP mixture (10 mM)	1.5 μ L
DTT (100 mM)	1.25 μ L
SuperScript III	1 μ L

- Incubate in a thermocycler:

- 25°C for 5 min
- 50°C for 1 hr
- 70°C for 15 min
- 4°C forever

- Purify cDNA with 45ul Agencourt AMPure XP beads. Elute with 25 μ L water. Transfer 23 μ L to a new tube.

3. First PCR

- Add the following to the 23 μ L cDNA from the previous step:

NEBNext (Q5) 2X master mix	25 μ L
Deg-PCR-1L primer (10 μ M)	1 μ L
Deg-PCR-1R primer (10 μ M)	1 μ L

- Run the following program in a PCR machine with a heated-lid:

1. 98°C for 30 sec
2. 98°C for 10 sec
3. 63°C for 10 sec
4. 72°C for 12 sec
5. Go to 2., 4 \times cycles total
6. 98°C for 10 sec
7. 70°C for 30 sec
8. 72°C for 12 sec
9. Go to 6., 6 \times cycles total
10. 72°C for 3 min
11. 11. 4°C forever

- Purify PCR products on 1% agarose gel. Run until 200- and 400bp ladder is separated well. Generally takes 45-60 mins when run at 100-120V. Use 100 bp DNA ladder to precisely select for 200–400 bp PCR products.
- It is important to cut just above 200nt ladder. To improve the visualization gel can be stained with EtBr wash.
- Run different genotypes on different gels to prevent cross contamination.
 - *It's normal that the PCR products are not visible on the gel at this point. Do not panic.*
- Extract PCR products with Qiagen gel extraction kit. Elute DNA with 24 μ L water. Transfer 23 μ L DNA to a new PCR tube.

4. Second PCR

- Add the following to the 23 μL DNA from the last step:

NEBNext (Q5) 2X master mix	25 μL
BPCR1 primer (10 μM)	1 μL
BPCRIdX primer (10 μM)	1 μL

- Run the following program in a PCR machine with a heated-lid:
 1. 98°C for 30 sec
 2. 98°C for 10 sec
 3. 70°C for 30 sec
 4. 72°C for 14 sec
 5. Go to 2., 4× cycles total
 6. 98°C for 10 sec
 7. 72°C for 34 sec
 8. Go to 6., 8× cycles total
 9. 72°C for 3 min
 10. 4°C forever
- Purify the final library with 50 μL Ampure XP beads. Elute with 27 μL EB and transfer 25 μL supernatant to a new tube.

5. Purification with non-denaturing PAGE:

- Prepare a 10% non-denaturing polyacrylamide gel (20 cm \times 16 cm \times 1.5 mm, length \times width \times thickness), with 20 small wells (each well: 17 mm \times 5 mm \times 1.5 mm, depth \times width \times thickness). The gel recipe is as follows:

water	32 mL
5x TBE	5 mL
40% Acrylamide/Bis-acrylamide (Sigma-Aldrich # A7802)	12.5
10% ammonium persulfate	500 μL
TEMED	50 μL

- Pre-run the gel in 0.5× TBE buffer at the constant power of 35 W for 30 min (the surface temperature of the gel reaches 45–55°C).
- To 25 μL of the library from the previous step add 5 μL of 6x Loading Dye (e.g., Trisack, Thermofisher #R1161) and mix thoroughly.
- Wash wells with the running buffer. Load 30 μL of libraries into wells keeping at least one well empty between samples. Remember to load 100bp DNA Ladder (e.g., NEB #N3231S) into the rightmost and the leftmost wells.
- Run the gel in 0.5× TBE at the constant power of 20 W for 75 min or until the xylene cyanol bands are just above the bottom of the gel.
- Open the glass plates and cut off one bottom corner of the gel to mark the loading sequence. Gently peel the gel off the glass plate and transfer it into the glass tray containing 1× SYBR Gold in 100 ml 0.5× TBE buffer. Stain the gel for 30 sec.
- Cut off margins of an Avery diamond clear standard weight sheet protector, separate the two pieces, and lay one piece on the bench with the inside facing up. Gently transfer the gel on the sheet protector. Always keep the gel on the sheet protector to avoid contamination.
- Scan the gel with Typhoon FLA 9000 or similar scanner.
- Transfer the gel with sheet protector on the UV box. Switch on the long-wave UV light and cut the gel between the **250bp** and **500bp** with a new razor blade (REMEMBER TO PROTECT YOUR EYES FROM UV LIGHT!!!). Do not use the

same blade to cut different samples. Transfer the gel slice into a 2.0 ml tube and crush it with either a single-use plastic pestel or by centrifuging it through a pierced 0.5ml tube as described in the Illumina Truseq Small RNA protocol. Add 0.7 ml of 0.3 M NaCl to each tube and rotate in cold room overnight to elute libraries.

- The following day, use a P1000 tip with the end cut off to transfer most of the gel elution mix from each 2.0 ml tube to a Spin-X column (Sigma #CLS8162), and spin at $10,000 \times g$ for 5 min to filter out gel pieces.
- Split elution from each 1.7 ml tube into 2 fresh tubes, $\sim 350 \mu\text{l}$ each, add $1 \mu\text{l}$ (15 μg) of GlycoBlue (Thermofisher # AM9515) to each tube and vortex briefly to mix.
- Add 3 volumes (1050 μl) of 100% ethanol to each tube and vortex briefly to mix. Incubate on ice for 1 hr.
- Spin at $20,000 \times g$ for 30 min at 4°C .
- Wash the pellet with 900 μl of 75% ethanol. Spin at $20,000 \times g$ for 5 min at 4°C and gently remove supernatant without touching the pellet.
- Spin at $20,000 \times g$ for 1 min at 4°C and gently remove the residual ethanol with a P10 tip.
- Air dry for 2 min and dissolve the pellet in 5–10 μl water in each tube and combine the library from the two tubes.
- Quantify libraries with Kapa

Recipes (Morlan et al. *PLOS ONE* 2012 & Adiconis et al. *Nature Methods* 2013)

rRNA oligo Hybridization buffer

100 mM Tris-Cl (pH 7.4)
200 mM NaCl

10X RNase H digestion buffer

500 mM Tris-Cl (pH 7.4)
1M NaCl
200 mM MgCl_2

10 \times 5' Ligation Buffer: 500 mM Tris-HCl (pH 7.8), 100 mM MgCl_2 , 100 mM DTT, and 10 mM ATP. Use freshly prepared buffer every prep as DTT and ATP degrade over time.

Water	10 μl
1 M Tris-HCl (pH 7.8)	50 μl
500 mM MgCl_2	20 μl
1 M DTT	10 μl
100mM rATP	10 μl
Final volume	100 μl